



zapA Antibody

Product Code	CSB-PA606126ZA01EYY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q11137
Immunogen	Recombinant Proteus mirabilis zapA protein
Raised In	Rabbit
Species Reactivity	Proteus mirabilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serralysin (EC 3.4.24.40) (Extracellular metalloprotease), zapA
Immunogen Species	Proteus mirabilis
Lead Time	Made-to-order (12-14 weeks)
Target Names	zapA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.