



TM1 Antibody

Product Code	CSB-PA375499ZA01ETF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1KYZ2
Immunogen	Recombinant Penaeus monodon TM1 protein
Raised In	Rabbit
Species Reactivity	Penaeus monodon
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tropomyosin (Tropomyosin, fast isoform) (Tm-Penm-fast), TM1
Immunogen Species	Penaeus monodon (Giant tiger prawn)
Lead Time	Made-to-order (12-14 weeks)
Target Names	TM1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.