



cyp c 1.02 Antibody

Product Code	CSB-PA20339ZA01EQE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8UUS2
Immunogen	Recombinant Cyprinus carpio cyp c 1.02 protein
Raised In	Rabbit
Species Reactivity	Cyprinus carpio
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Parvalbumin (cyp c 1.02)
Immunogen Species	Cyprinus carpio (Common carp)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	cyp c 1.02
Product Modify	Unmodified