



# GOT2 Antibody

<b>Product Code</b>	CSB-PA009681ZA01BO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P12344
<b>Immunogen</b>	Recombinant Bos taurus GOT2 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bos taurus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein A/G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Aspartate aminotransferase, mitochondrial (mAspAT) (EC 2.6.1.1) (Fatty acid-binding protein) (FABP-1) (Glutamate oxaloacetate transaminase 2) (Kynurenine aminotransferase 4) (Kynurenine aminotransferase IV) (Kynurenine--oxoglutarate transaminase 4) (Kynurenine--oxoglutarate transaminase IV) (Plasma membrane-associated fatty acid-binding protein) (FABPpm) (Transaminase A), GOT2
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	GOT2
<b>Product Modify</b>	Unmodified