



CDC48 Antibody

Product Code	CSB-PA346583ZA01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54774
Immunogen	Recombinant Glycine max CDC48 protein
Raised In	Rabbit
Species Reactivity	Glycine max
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division cycle protein 48 homolog (Valosin-containing protein homolog) (VCP), CDC48
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CDC48
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.