



Ion Antibody

Product Code	CSB-PA305230ZA01MLW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78025
Immunogen	Recombinant Mycoplasma pneumoniae Ion protein
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lon protease (EC 3.4.21.53) (ATP-dependent protease La), lon
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	lon
Product Modify	Unmodified