



Glp1r Antibody

Product Code	CSB-PA009514ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32301
Immunogen	Recombinant Rattus norvegicus Glp1r protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Glucagon-like peptide 1 receptor (GLP-1 receptor) (GLP-1-R) (GLP-1R), Glp1r, Glpr
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Glp1r
Product Modify	Unmodified