



bsk Antibody

Product Code	CSB-PA310094ZA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P92208
Immunogen	Recombinant Drosophila melanogaster bsk protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stress-activated protein kinase JNK (dJNK) (EC 2.7.11.24) (Protein basket), bsk, JNK
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (12-14 weeks)
Target Names	bsk
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.