



Glutenin, high molecular weight subunit PC237 Antibody

Product Code	CSB-PA355887ZA01TQN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02862
Immunogen	Recombinant Triticum aestivum Glutenin, high molecular weight subunit PC237 protein
Raised In	Rabbit
Species Reactivity	Triticum aestivum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutenin, high molecular weight subunit PC237 (Fragment)
Immunogen Species	Triticum aestivum (Wheat)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified