



ORF3 Antibody

Product Code	CSB-PA325916ZA01PQU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17521
Immunogen	Recombinant Potato leafroll virus ORF3 protein
Raised In	Rabbit
Species Reactivity	Potato leafroll virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major capsid protein (Coat protein) (CP)
Immunogen Species	Potato leafroll virus (strain Potato/Canada/Rowhani/1979) (PLrV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ORF3
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.