



RHOT1 Antibody

Product Code	CSB-PA720036ZA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ZM73
Immunogen	Recombinant Gallus gallus RHOT1 protein
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial Rho GTPase 1 (MIRO-1) (EC 3.6.5.-) (Ras homolog gene family member T1), RHOT1
Immunogen Species	Gallus gallus (Chicken)
Lead Time	Made-to-order (12-14 weeks)
Target Names	RHOT1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.