



DMP1 Antibody

Product Code	CSB-PA20439ZA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F1PHT7
Immunogen	Recombinant Canis lupus familiaris DMP1 protein
Raised In	Rabbit
Species Reactivity	Canis lupus familiaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dentin matrix acidic phosphoprotein 1 (DMP1)
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	DMP1
Product Modify	Unmodified