



PRDX1 Antibody

Product Code	CSB-PA872876ZA01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JKY1
Immunogen	Recombinant Cricetulus griseus PRDX1 protein
Raised In	Rabbit
Species Reactivity	Cricetulus griseus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxiredoxin-1 (EC 1.11.1.15) (Thioredoxin peroxidase 2) (TPX-2), PRDX1, TDPX2
Immunogen Species	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PRDX1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.