



h1-0-a Antibody

Product Code	CSB-PA332862ZA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22845
Immunogen	Recombinant Xenopus laevis h1-0-a protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H1.0-A (H1-SB) (H1E) (Histone H1(0)-1) (Histone H5B) (XIH5B), h1f0-a
Immunogen Species	Xenopus laevis (African clawed frog)
Lead Time	Made-to-order (12-14 weeks)
Target Names	h1-0-a
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.