



pepA Antibody

Product Code	CSB-PA302719ZA01MLW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P75206
Immunogen	Recombinant Mycoplasma pneumoniae pepA protein
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable cytosol aminopeptidase (EC 3.4.11.1) (Leucine aminopeptidase) (LAP) (EC 3.4.11.10) (Leucyl aminopeptidase), pepA
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pepA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.