



Dpp4 Antibody

Product Code	CSB-PA007139ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14740
Immunogen	Recombinant Rattus norvegicus Dpp4 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidyl peptidase 4 (EC 3.4.14.5) (Bile canaliculus domain-specific membrane glycoprotein) (Dipeptidyl peptidase IV) (DPP IV) (GP110 glycoprotein) (T-cell activation antigen CD26) (CD antigen CD26) [Cleaved into: Dipeptidyl peptidase 4 membrane form (Dipeptidyl peptidase IV membrane form); Dipeptidyl peptidase 4 soluble form (Dipeptidyl peptidase IV soluble form); Dipeptidyl peptidase 4 60 kDa soluble form (Dipeptidyl peptidase IV 60 kDa soluble form)], Dpp4, Cd26
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Dpp4
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.