



## E Antibody

Product Code	CSB-PA360435ZA01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00720
Immunogen	Recombinant Enterobacteria phage T4 E protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endolysin (EC 3.2.1.17) (Lysis protein) (Lysozyme) (Muramidase), E
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	E
Product Modify	Unmodified