



NP Antibody

Product Code	CSB-PA630945ZA01IFY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1PUD5
Immunogen	Recombinant Influenza A virus NP protein
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), NP
Immunogen Species	Influenza A virus (strain A/Port Chalmers/1/1973 H3N2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	NP
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.