



groEL Antibody

Product Code	CSB-PA314770ZA01DSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0Z7
Immunogen	Recombinant Chlamydia trachomatis groEL protein
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	60 kDa chaperonin (57 kDa chlamydial hypersensitivity antigen) (GroEL protein) (Heat shock protein 60) (HSP60) (Protein Cpn60), groL, groEL hypB mopA
Immunogen Species	Chlamydia trachomatis (strain D/UW-3/Cx)
Lead Time	Made-to-order (12-14 weeks)
Target Names	groEL
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.