



# Becn1 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA002658ZA01MO  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | O88597  |
| <b>Immunogen</b>           | Recombinant Mus musculus Becn1 protein  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Mus musculus  |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4  |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Beclin-1 (Coiled-coil myosin-like BCL2-interacting protein) [Cleaved into: Beclin-1-C 35 kDa; Beclin-1-C 37 kDa], Becn1 |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)  |
| <b>Lead Time</b>           | Made-to-order (12-14 weeks)   |
| <b>Gene Names</b>          | Becn1   |
| <b>Product Modify</b>      | Unmodified  |