




wac Antibody

Product Code	CSB-PA319157ZA01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10104
Immunogen	Recombinant Enterobacteria phage T4 Fibrin (wac) (2-487AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Fibrin (Collar protein) (Whisker antigen control protein), wac
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Target Names	wac

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: wac antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67.7 kDa

Observed band size: 68 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.