



hspX Antibody

Product Code	CSB-PA358651ZA01MVH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A5B8
Immunogen	Recombinant Mycobacterium bovis hspX protein
Raised In	Rabbit
Species Reactivity	Mycobacterium bovis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	14 kDa antigen (16 kDa antigen) (HSP 16.3), hspX
Immunogen Species	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
Lead Time	Made-to-order (12-14 weeks)
Target Names	hspX
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.