



# isdB Antibody

Product Code	CSB-PA313758ZA01SUM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C7J5
Immunogen	Recombinant Staphylococcus aureus isdB protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-regulated surface determinant protein B (Fur-regulated protein B) (Staphylococcal iron-regulated protein H) (Staphylococcus aureus surface protein J), isdB, frpB sasJ sirH
Immunogen Species	Staphylococcus aureus (strain USA300 / TCH1516)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	isdB
Product Modify	Unmodified