



gag Antibody

Product Code	CSB-PA724163ZA01GBAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q67652
Immunogen	Recombinant Giardia lamblia virus gag protein
Raised In	Rabbit
Species Reactivity	Giardia lamblia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major capsid protein (Gag protein) (Major coat protein), gag
Immunogen Species	Giardia lamblia virus (isolate Wang) (GLV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	gag
Product Modify	Unmodified