



IZUMO1R Antibody

Product Code	CSB-PA008789ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6ND01
Immunogen	Recombinant Homo sapiens IZUMO1R protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm-egg fusion protein Juno (Folate receptor 4) (Folate receptor delta) (FR-delta) (IZUMO1 receptor protein JUNO), IZUMO1R, FOLR4 JUNO
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IZUMO1R
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.