



eaeA Antibody

Product Code	CSB-PA346794ZA01HAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52869
Immunogen	Recombinant Hafnia alvei eaeA protein
Raised In	Rabbit
Species Reactivity	Hafnia alvei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intimin (Attaching and effacing protein) (Outer membrane protein) (Fragment), eaeA, eae
Immunogen Species	Hafnia alvei
Lead Time	Made-to-order (12-14 weeks)
Gene Names	eaeA
Product Modify	Unmodified