



PNC1 Antibody

Product Code	CSB-PA346834ZA01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53184
Immunogen	Recombinant Saccharomyces cerevisiae PNC1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotinamidase (EC 3.5.1.19) (Nicotinamide deamidase) (NAMase), PNC1
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	PNC1
Product Modify	Unmodified