



glpE Antibody

Product Code	CSB-PA533026ZA01ENP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1IP41
Immunogen	Recombinant Escherichia coli glpE protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Thiosulfate sulfurtransferase GlpE (EC 2.8.1.1), glpE
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	glpE
Product Modify	Unmodified