



yscM Antibody

Product Code	CSB-PA300275ZA01YAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69979
Immunogen	Recombinant Yersinia pseudotuberculosis serotype I yscM protein
Raised In	Rabbit
Species Reactivity	Yersinia pseudotuberculosis serotype I
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Yop proteins translocation protein M (Low calcium response locus protein Q), yscM, lcrQ
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Lead Time	Made-to-order (12-14 weeks)
Target Names	yscM
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.