



CMA1 Antibody

Product Code	CSB-PA005599ZA01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56435
Immunogen	Recombinant Macaca fascicularis CMA1 protein
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chymase (EC 3.4.21.39) (Alpha-chymase), CMA1
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	CMA1
Product Modify	Unmodified