



catA Antibody

Product Code	CSB-PA357202ZA01AWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07773
Immunogen	Recombinant Acinetobacter baylyi catA protein
Raised In	Rabbit
Species Reactivity	Acinetobacter baylyi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Catechol 1,2-dioxygenase (EC 1.13.11.1) (1,2-CTD), catA
Immunogen Species	Acinetobacter baylyi
Lead Time	Made-to-order (12-14 weeks)
Target Names	catA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.