



AER Antibody

Product Code	CSB-PA654400ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q39172
Immunogen	Recombinant Arabidopsis thaliana AER protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADPH-dependent oxidoreductase 2-alkenal reductase (AtAER) (EC 1.3.1.-) (EC 1.3.1.74) (NADP-dependent alkenal double bond reductase P1) (DBR1) (NADPH-azodicarbonyl/quinone reductase) (NADPH:2-alkenal/one alpha,beta-hydrogenase) (ALH) (P1-zeta-crystallin protein) (P1-ZCr), AER, P1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (12-14 weeks)
Target Names	AER
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.