



dps Antibody

Product Code	CSB-PA337381ZA01HUV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43313
Immunogen	Recombinant Helicobacter pylori dps protein
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA protection during starvation protein (EC 1.16.-.-) (Bacterioferritin) (HP-NAP) (Neutrophil-activating protein A) (NAP A), dps, napA
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Lead Time	Made-to-order (12-14 weeks)
Target Names	dps
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.