



Ipp Antibody

Product Code	CSB-PA355906ZA01EMS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02939
Immunogen	Recombinant Erwinia amylovora Ipp protein
Raised In	Rabbit
Species Reactivity	Erwinia amylovora
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane lipoprotein (Murein-lipoprotein), Ipp
Immunogen Species	Erwinia amylovora
Lead Time	Made-to-order (12-14 weeks)
Target Names	Ipp
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.