



U1-cyrtautoxin-As1c Antibody

Product Code	CSB-PA344579ZA01AMR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49270
Immunogen	Recombinant Apomastus schlingeri U1-cyrtautoxin-As1c protein
Raised In	Rabbit
Species Reactivity	Apomastus schlingeri
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	U1-cyrtautoxin-As1c (U1-CUTX-As1c) (Aptotoxin VI) (Aptotoxin-6) (Paralytic peptide VI) (PP VI)
Immunogen Species	Apomastus schlingeri
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.