



omcA Antibody

Product Code	CSB-PA316432ZA01DSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CC05
Immunogen	Recombinant Chlamydia trachomatis omcA protein
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small cysteine-rich outer membrane protein OmcA (Small-CRP) (9 kDa cysteine-rich lipoprotein) (9kDa-CRP), omcA, omp2A omp3
Immunogen Species	Chlamydia trachomatis (strain D/UW-3/Cx)
Lead Time	Made-to-order (12-14 weeks)
Target Names	omcA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.