



Aspergillopepsin-2 Antibody

Product Code	CSB-PA338616ZA01DOZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24665
Immunogen	Recombinant Aspergillus niger Aspergillopepsin-2 protein
Raised In	Rabbit
Species Reactivity	Aspergillus niger
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aspergillopepsin-2 (EC 3.4.23.19) (Acid protease A) (Aspergillopepsin II) (Proctase A) [Cleaved into: Aspergillopepsin-2 light chain (Aspergillopepsin II light chain); Aspergillopepsin-2 heavy chain (Aspergillopepsin II heavy chain)]
Immunogen Species	Aspergillus niger
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.