



LHB Antibody

Product Code	CSB-PA012910ZA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18842
Immunogen	Recombinant Canis lupus familiaris LHB protein
Raised In	Rabbit
Species Reactivity	Canis lupus familiaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lutropin subunit beta (Luteinizing hormone subunit beta) (LH-B) (LSH-B) (LSH-beta) (Lutropin beta chain) (Fragment), LHB
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	LHB
Product Modify	Unmodified