



MSRB7 Antibody

Product Code	CSB-PA848741ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VY86
Immunogen	Recombinant Arabidopsis thaliana MSRB7 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptide methionine sulfoxide reductase B7 (AtMSRB7) (EC 1.8.4.12) (Peptide-methionine (R)-S-oxide reductase), MSRB7
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (12-14 weeks)
Target Names	MSRB7
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.