



Lgmn Antibody

Product Code	CSB-PA012903ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O89017
Immunogen	Recombinant Mus musculus Lgmn protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Legumain (EC 3.4.22.34) (Asparaginyl endopeptidase) (Protease, cysteine 1), Lgmn, Prsc1
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Lgmn
Product Modify	Unmodified