



UL83 Antibody

Product Code	CSB-PA361930ZA01HWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06725
Immunogen	Recombinant Human cytomegalovirus UL83 protein
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	65 kDa phosphoprotein (pp65) (65 kDa matrix phosphoprotein) (Phosphoprotein UL83) (Tegument protein UL83), UL83
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (Human herpesvirus 5)
Lead Time	Made-to-order (12-14 weeks)
Target Names	UL83
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.