



peb1A Antibody

Product Code	CSB-PA379898ZA01CWF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1VZQ4
Immunogen	Recombinant Campylobacter jejuni subsp. jejuni serotype O:23/36 peb1A protein
Raised In	Rabbit
Species Reactivity	Campylobacter jejuni subsp. jejuni serotype O:23/36
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major cell-binding factor (CBF1) (PEB1), peb1A, pebA
Immunogen Species	Campylobacter jejuni subsp. jejuni serotype O:23/36 (strain 81-176)
Lead Time	Made-to-order (12-14 weeks)
Target Names	peb1A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.