



groEL Antibody

Product Code	CSB-PA656857ZA01SGF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q328C4
Immunogen	Recombinant Shigella dysenteriae serotype 1 groEL protein
Raised In	Rabbit
Species Reactivity	Shigella dysenteriae serotype 1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	60 kDa chaperonin (GroEL protein) (Protein Cpn60), groL, groEL
Immunogen Species	Shigella dysenteriae serotype 1 (strain Sd197)
Lead Time	Made-to-order (12-14 weeks)
Target Names	groEL
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.