



ZNF160 Antibody

Product Code	CSB-PA026549ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCG1
Immunogen	Recombinant Homo sapiens ZNF160 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 160 (Zinc finger protein HZF5) (Zinc finger protein Kr18) (HKr18), ZNF160, KIAA1611
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	ZNF160
Product Modify	Unmodified