




fimA Antibody

Product Code	CSB-PA361210ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04128
Immunogen	Recombinant Escherichia coli Type-1 fimbrial protein, A chain (fimA) (24-182AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Type-1 fimbrial protein, A chain (Type-1A pilin), fimA, pilA
Immunogen Species	Escherichia coli (strain K12)
Target Names	fimA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: fimA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31.8 kDa

Observed band size: 36 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.