



FGF10 Antibody

Product Code	CSB-PA008616ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15520
Immunogen	Recombinant Human Fibroblast growth factor 10 (FGF10), partial (37-208AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 10 (FGF-10) (Keratinocyte growth factor 2), FGF10
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Target Names	FGF10
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.