



emm5 Antibody

Product Code	CSB-PA355917ZA01SMY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02977
Immunogen	Recombinant Streptococcus pyogenes serotype M5 emm5 protein
Raised In	Rabbit
Species Reactivity	Streptococcus pyogenes serotype M5
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	M protein, serotype 5, emm5, smp5
Immunogen Species	Streptococcus pyogenes serotype M5
Lead Time	Made-to-order (12-14 weeks)
Target Names	emm5
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.