



cdtB Antibody

Product Code	CSB-PA671598ZA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q46669
Immunogen	Recombinant Escherichia coli cdtB protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytolethal distending toxin subunit B (CDT B) (Deoxyribonuclease CdtB) (EC 3.1.-.-), cdtB
Immunogen Species	Escherichia coli
Lead Time	Made-to-order (12-14 weeks)
Target Names	cdtB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.