



CYND1 Antibody

Product Code	CSB-PA516895ZA01EQB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O04701
Immunogen	Recombinant Cynodon dactylon CYND1 protein
Raised In	Rabbit
Species Reactivity	Cynodon dactylon
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Cyn d 1 (allergen Cyn d 1), CYND1
Immunogen Species	Cynodon dactylon (Bermuda grass) (Panicum dactylon)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CYND1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.