



Sult1a1 Antibody

Product Code	CSB-PA022933ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17988
Immunogen	Recombinant Rattus norvegicus Sult1a1 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfotransferase 1A1 (ST1A1) (EC 2.8.2.1) (Aryl sulfotransferase) (Aryl sulfotransferase IV) (ASTIV) (Minoxidil sulfotransferase) (Mx-ST) (PST-1) (Phenol sulfotransferase) (Sulfokinase) (Tyrosine-ester sulfotransferase), Sult1a1, St1a1
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Sult1a1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.